



**FILTER ASSEMBLY,  
MEDICAL LAB APPLICATION**

**SANITARY FABRICATION AND FINISH**

**HAND DISASSEMBLY/REASSEMBLY**

**EASY CLEANING**

**Features:**

- Cover, cone with view port
- Gasket to suite
- Baffled mixing chamber
- Screen (sieve) welded assembly with mesh to suite
- Drain piped to outside container
- Base attachment to vibration table
- Materials: SS316
- Finish, inside & outside, 240 grit

**NOTE: TOP ENTRY MIXER IS NOT  
SHOWN**



**DANTCO MIXERS CORPORATION  
9 OAK STREET**

**PATERSON, NJ 07501-3013**

**973-278-8776, FAX 973-278-1538**

**EMAIL [sales@dantco.com](mailto:sales@dantco.com), WEB <http://www.dantco.com>**